EDENTA AG, Hauptstrasse 7, CH-9434 AU/SG – **Switzerland**
Tel.: +41 71 747 25 25 Fax: +41 71 747 25 50 e-mail: info@edenta.ch

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GB User Information for rotating dental and surgical instruments

- Use of improper instruments can cause insufficient preparation results.
- Excessive speeds while using instruments with overly large heads endanger the stability of the instruments.
- Excessive speeds applied near to the dental pulp can endanger its vitality. The same can be the result when using an inadequate supply of water spray.
- Observe permissible speeds indicated on packaging.
- The handpiece must be in proper working condition. The use of an adequate supply of water spray of a minimum of 50 ml/min on the working part is necessary in the 1500 min⁻¹ speed range to ensure reliable cooling of the material being worked on.
- Additional external cooling is needed for instruments with a total length of > 19 mm and a head diameter of 1.8 mm (ISO -018).
- Moistening of the surfaces to be prepared must not be impaired by an incorrect suction system. Deflection of the coolant jet can also cause heat related damage.

Detailed recommendations for preparation in accordance with DIN EN ISO 17664

www.edenta.com

- Sterilisation for polishers at 134°C / 2.1 bar; 5 min. is recommended.
 - The use of bur stands or trays is necessary to avoid contact damage.
 - The use of inappropriate disinfectants and/or keeping the instruments too long in the disinfectant can result in corrosion (increased danger of breakage) and discolouration of the instruments or the colour coding can come off. No alkaline disinfectant must be used for the disinfection of silicon rubber polishers.
 - Do not place polishers in alcohol!
 - Use brushes / Flexi-Snap as single-use products! ☒
- Steel instruments** are unsuitable for both mechanical preparation and the steam steriliser. They should be replaced by corresponding tungsten-carbide instruments. **Single-use products** (identified on the packaging with ☒) are not approved for reuse. Safe use cannot be guaranteed if these products are reused, as there is the risk of infection and/or the products are no longer safe to use. The end of the product's service life is generally determined by wear and damage due to use.

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- The recommended applied contact pressure range is 0.3 – 2 N (30 – 200 p) depending on the kind of preparation.
- Avoid tilting, levering, or blocking the instruments up by using too much contact pressure (increased danger of fracture).
- Insert the instruments as deeply as possible without excessive use of force into the chuck and check to ensure that they are secure.
- Extremely tapered instruments will wear faster at the tip due to the smaller size. To maintain the longevity of these instruments, careful preparation with less contact pressure is important.
- Apply Turbo instruments over their entire grinding surface.
- Roughness obtained with a diamond coating >ISO 524 has to be finished.
- Immediately remove any instruments that are damaged, bent or no longer run concentrically.
- Depending on the application, the use of protective goggles is recommended.
- Rotating instruments have to be disinfected, cleaned, dried and sterilised before the first application on the patient and immediately after every use.

Storage and Keeping Conditions:

In dry conditions and protected against contaminants. Protect instruments in general against chemicals, acids, heat and extreme temperature variations. Do not use polishers after expiration date indicated.

Observe guidelines from Robert-Koch-Institut - ch. E3, (Hygiene Requirements in Dentistry) or the valid legal and hygiene regulations in the respective country!

Steel instruments are unsuitable for both mechanical preparation and the steam steriliser. They should be replaced by corresponding tungsten-carbide instruments. **Single-use products** (identified on the packaging with ☒) are not approved for reuse. Safe use cannot be guaranteed if these products are reused, as there is the risk of infection and/or the products are no longer safe to use. The end of the product's service life is generally determined by wear and damage due to use.

Colour-anodised aluminium units (e.g. Bur block 40500 to 40580 and Retopin mandrel) lose their colour during use of standard cleaning procedures and in a cleaning and disinfection machine. Cleaning and disinfectant agents, which are specifically designed for these materials should be used during preparation.

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