EDENTA AG, Hauptstrasse 7, CH-9434 AU/SG – Switzerland
Tel.: +41 71 747 25 25 Fax: +41 71 747 25 50 e-mail: info@edenta.ch

EDENTA Etablissement, Industriestrasse 13, LI-9486 Schaanwald – Liechtenstein
Tel.: +423 375 20 50 Fax: +423 375 20 55 e-mail: info@edenta.com

DE

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GB

Product Information and Instructions for Diamond Finishing Strips:

Description of the product:

Diamond finishing strips are metal strips diamond-coated on one side. They are available in three different grits and used chairside for breaking the interproximal contacts.

Range of indications:

- Trimming excess material on the proximal regions of fillings.
- Contouring the proximal regions of all types of filling.
- For use on composites, acrylics, glass-ionomer cements, amalgams etc.
- Removing excess cement after placing inlays, onlays and laminate veneers.

Application:

Non-serrated diamond finishing strips:

As the center section of these strips is not diamond-coated, they can be pulled backwards and forwards through embrasures. The diamond-coated surface must point toward the filling being trimmed. The diamond-coated section is pulled over the filling to trim off any excess material and contour the filling anatomically. The strips are used in descending order, from medium grit down to superfine.

They are colour-coded to indicate their grit sizes:

Blue	Medium grit
Red	Fine grit
Yellow	Superfine grit

Important, please note:

When pulling the strip backwards and forwards, ensure that its edge does not injure the interdental papillae. This can be prevented by tilting the strip to contact the tooth closely after sliding it into the embrasure. Never kink the strip as this would damage it and split the diamond-coating off.

Serrated diamond finishing strips:

The serrated edge is used for separating the interdental spaces. Once the serrated edge has been slipped through the contact areas, turn the strip over and continue as described above. This avoids injury to the papillae.

Disinfection, cleaning and sterilization:


Diamond finishing strips must always be disinfected, cleaned, dried and sterilized prior to first use and after use in patients' mouths.

Detailed recommendations for preparation in accordance with DIN EN ISO 17664

www.edenta.com

Disinfection:

Never kink the strip while disinfecting, cleaning and sterilizing. Only use disinfectants approved by the German Society for Hygiene and Microbiology for use on rotary instruments. Adhere strictly to

the disinfectant manufacturer's instructions. Once the prescribed reaction time has elapsed, remove the diamond finishing strip from the disinfectant and, if prescribed, rinse with running water and dry.

Cleaning:

If the diamond strip is cleaned in an ultrasonic unit, never place it in the same vessel as other instruments, especially diamond instruments, as this may damage the strip. The diamond strips can also be cleaned with a brush - avoid damaging the diamond-coating.

Sterilization:

Sterilize the diamond strips in an autoclave as usual at 134°C and 2.1 bars for five minutes.

Checking:

Before placing these strips in storage, always check them for damage. If the diamond-coating is worn or damaged, or the strip kinked, discard it.

Storage:

Keep in a dry place, protected against recontamination. Do not kink or damage in any other manner.

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