EDENTA AG, Hauptstrasse 7, CH-9434 AU/SG – Switzerland

Tel.: +41 71 747 25 25 Fax: +41 71 747 25 50 e-mail: info@edenta.ch

EDENTA Etabl., Industriest. 13, LI-9486 Schaanwald – Liechtenstein

Tel.: +423 375 20 50 Fax: +423 375 20 55 e-mail: info@edenta.com

DE

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GB

Product information and instructions for use for podiatry instruments:

Scope:

Rotary instruments for medical (podiatry) and cosmetic footcare and for nail design.

Uses:

Medical footcare (podiatry)

- Treatment of curling and ingrowing nails from fungal nail infections or thickened nails.
- Removal of excess hard skin and calluses.
- Correct removal and treatment of corns and warts.

Cosmetic footcare

- Correctly carrying out skincare and decorative procedures on healthy feet.

Nail design

- Treatment of artificial nails.

Use:

All podiatry instruments were developed and constructed for specific use. Incorrect use can therefore lead to damage to tissues, premature wear, destruction of the instruments and to risk for the user, patients and third parties.

- It should be ensured that only technically and hygienically faultless, maintained and cleaned turbines and handpieces are used.
- The instruments must be inserted as deeply as possible.
- Podiatry instruments should reach the correct speed before being applied to the object.
- Tilting or levering must be avoided; this leads to an increased risk of breakage.
- Use a suction or spray device during milling, grinding and polishing on account of the dust production. Wearing protective glasses and facemask is recommended.
- Use disposable gloves when handling the instruments as otherwise there is a risk of infection and patient safety is no longer met.
- To avoid undesirable heat, use a low speed or adequate water cooling during treatment.
- Excessive pressure can lead to damage to skin and tissue due to overheating.
- When the diameter of the working part is greater than the shaft thickness, strong centrifugal forces can occur at excessive speeds, which can lead to bending of the shaft and / or breakage of the instruments. For this reason the maximum permitted speed must never be exceeded. You will find the maximum permitted speed in the information on the packaging.

Disinfection, cleaning and sterilisation:

Before using on a patient for the first time and after every use the instruments must be disinfected, cleaned, dried and, where necessary, sterilised.

Detailed recommendations for preparation in accordance with DIN EN ISO 17664

Storage / transport:

Immediately after use, submerge the instruments completely in a suitable (alkali- and aldehyde-free) cleaning/disinfectant solution, e.g., HELVEMED Instrument Forte and transfer to reprocessing after one hour at the latest. Follow the manufacturer's instructions regarding concentration and exposure time.

Processing procedure:

1. Remove the instruments from the disinfectant bath and rinse off surface soiling thoroughly with clean running water. Rinse off adherent soiling of the instruments with a hard plastic brush under clean running water. Note: use of an ultrasonic device is recommended for cleaning instruments with interstices. When cleaning in an ultrasonic bath, the instruments must not come in contact with each other so they should be placed in a suitable instrument stand in the ultrasonic device filled with cleaning and disinfectant agents. Since the vibrations in the ultrasonic bath can be absorbed by the materials of the polishers, these should be processed only in an instrument bath. For cleaning and chemical disinfection in the ultrasonic device, follow the cleaning/disinfectant agent manufacturer's instructions regarding concentration and exposure time. The exposure time starts only when the last instrument has been placed in the ultrasonic device and must not be reduced. After the exposure time is complete, rinse off any residual disinfectant on the instruments thoroughly with clean running water and dry them. Inspect visually for cleanliness (e.g., using a loupe etc. with 8 to 10-fold magnification). If residual contamination is identified, repeat cleaning and chemical disinfection until no further contamination is visible.
2. Inspect the cleaned and dried instruments for intactness (e.g., using a loupe etc. with 8 to 10-fold magnification) and reject if necessary (absence of diamond coating, blunt/broken cutting edges, fractured working parts, corroded surfaces, bent instruments etc.)
3. Terminal sterilisation at 134°C, 2.1 bar, 5 minutes, in an autoclave is recommended. To avoid staining and corrosion, the steam must be free from substances. The instruments should be sealed in suitable films prior to sterilisation. The device manufacturer's instructions should be followed.

Storage:

Storage of rotary instruments should be at room temperature in processed and hygienically maintained stands, bowls or other suitable containers and in the original packaging until first use. The same applies for sterilised instruments and those packed for sterilisation. They must be stored protected from dust, moisture and recontamination.

Material stability:

When selecting cleaning and disinfectant agents, ensure that they do not contain the following components: - organic, mineral and oxidising acids / - strong alkalis (pH > 10.5 not permissible, only neutral or weak alkaline cleaning agents are recommended) / - do not use alkaline cleaning agents for polishers / - alcohols, ethers and ketones, petroleum / - oxidising agents. Never clean instruments and sterilisation trays with metal brushes or steel wool. Colour-anodised aluminium components lose their colour when conventional cleaning methods are used. Special cleaning and disinfectant agents (e.g., HELVEMED Instrument Forte) should be used for processing.

Reprocessing limitation:

The end of the product's useful life is determined essentially by wear and damage due to use. Frequent reprocessing has no influence on the instruments' performance. Disposable items (labelled on the packaging) are not authorised for reuse. Non-hazardous use cannot be guaranteed if these products are used again since there is a risk of infection and/or the safety of the products is no longer assured.

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