EDENTA AG, Hauptstrasse 7, CH-9434 AU/SG – **Switzerland**
Tel.: +41 71 747 25 25 Fax: +41 71 747 25 50 e-mail: info@edenta.ch

EDENTA Etabl., Industriestr. 13, LI-9486 Schaanwald – **Liechtenstein**
Tel.: +423 375 20 50 Fax: +423 375 20 55 e-mail: info@edenta.com

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User Information for NT Endo Files (engine use)

Indication and efficacy

- The NT-Endo system was developed to meet the requirements for engine use root canal preparation according to the crown-down technique. The specially adapted blade geometry of the nickel-titanium instruments guarantees smooth, conically prepared root canals. Do not use this instrument for any purpose except those applications listed above.
- For use by dentists only.

Precautions for use

- To prevent infection, autoclave the product (134°C / 2.1 bar; 5 min.) and ensure that it is fully sterilized before use.
- Adhere strictly to the maximum permissible speed of 300 r.p.m. as exceeding this speed may cause fracturing.
- Use an endomotor-handpiece with controllable torque.
- Maximum torque for .04 taper version: 1.5 - 1.8 Ncm. Maximum torque for .06 taper version: 2.4 - 2.9 Ncm.
- Adhere strictly to the handpiece manufacturer's instructions and insert the shank into the chuck fully to prevent it jamming half way.
- Place a rubber dam etc. in the patient's mouth.
- Take an X-ray to confirm the working length.

- Choose the most appropriate type of instrument for each treatment and follow the generally accepted procedures of crown-down technique.
- Before using, examine the instruments outside of the oral cavity for deformation, scratches and cracks.
- Before using, run the instruments outside of the oral cavity to make sure that they are true-running.
- As products with thin, long or large working sections may fracture or twist, avoid using them at an inappropriate angle and exerting excessive pressure.

Maintenance and inspections

- Caution, as these instruments has sharp edges, it must be used with great care to avoid injury to fingers.
- After using, wash the instrument with medical cleaning agent and brush, then wash away foreign substances like adherent body fluids and body tissues completely and sterilize. When using medical cleaning agent, adhere strictly to its manufacturer's instruction manual.
- To avoid damaging its working parts, the instruments must be placed on a stand for cleaning in an ultrasonic cleaner.
- After cleaning, the instruments must be sterilized in an autoclave under the following conditions: 134°C / 2.1 bar; 5 min.
- Do not use a high pressure steam sterilizer which heats to more than 200°C, including during the drying process.
- Dispose of the instrument if damaged or its performance decreases.
- The instrument should be disposed of as if it were infectious medical waste.

Detailed recommendations for preparation in accordance with DIN EN ISO 17664 www.edenta.com

Storage and duration of use

- Avoid storing at high temperature, in humid conditions and direct sunlight. Keep away from liquids. Store at room temperature.
- Enforce the first-in first-out method for stock management.



Crown-Down-Shaping technique

Canal sizes

Three different sets of instruments are available to match the size of the canal; narrow, medium or wide.

Preparation sequence (crown-down technique)

1. Prepare straight coronal canal access region using file No. 10.
2. Create the guide path using a No. 010/015 K manual file and determine the working length (WL).
3. Then prepare to approx. 3 mm short of the working length (WL) using a No. 06 and No. 04 file.
4. Apical finishing then follows using a No. 06 and No. 04 file.

Small canal NT-901SO				
NT	10	06	04	06
Orifice shaping				
10% - 35				
Crown-Down				
6% - 25				
4% - 25				
6% - 20				
4% - 20				
Working Length determination				
K-Files 010/015				
Apical finishing				
6% - 20				
4% - 25				
Final flaring if necessary				
6% - 25				

Medium canal NT-902SO				
NT	10	06	04	06
Orifice shaping				
10% - 35				
Crown-Down				
6% - 30				
4% - 30				
6% - 25				
4% - 25				
Working Length determination				
K-Files 010/015				
Apical finishing				
6% - 25				
4% - 30				
Final flaring if necessary				
6% - 30				

Large canal NT-903SO				
NT	10	06	04	06
Orifice shaping				
10% - 35				
Crown-Down				
6% - 35				
4% - 35				
6% - 30				
4% - 30				
Working Length determination				
K-Files 010/015				
Apical finishing				
6% - 30				
4% - 35				
Final flaring if necessary				
6% - 35				

ST-Reamer 10%				
Taper:	.10			
Length:	21 mm			
Size:				35
Art.-No.:				ST-10.35

NT-Files 4%				
Taper:	.04			
Length:	25 mm			
Size:	20	25	30	35
Art.-No.:	NT-04.20	NT-04.25	NT-04.30	NT-04.35

NT-Files 6%				
Taper:	.06			
Length:	25 mm			
Size:	20	25	30	35
Art.-No.:	NT-06.20	NT-06.25	NT-06.30	NT-06.35

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