EDENTA AG, Hauptstrasse 7, CH-9434 AU/SG – Switzerland
Tel.: +41 71 747 25 25 Fax: +41 71 747 25 50 e-mail: info@edenta.ch

EDENTA Etabl., Industriestr. 13, LI-9486 Schaanwald – Liechtenstein
Tel.: +423 375 20 50 e-mail: info@edenta.com

DE Anwendungshinweise Head Master, active:

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GB Head Master, active Instruction for Use:

For Direct Restoration With Plastic Filling Materials. System has to be disinfected, cleaned, dried and sterilised before the first use; with steam (autoclave), temperature of 134°C, pressure of 2.1 bar, exposure time 5 min.

Detailed recommendations for preparation in accordance with DIN EN ISO 17664 www.edenta.com

- a) Treat the root of the tooth endodontically.
 - b) Remove as much the filling of the root as possible.
 - c) With the root canal drill prepare the root.
 - d) Drill the canal to the desired depth. Recommended speed: approx. 500-1000 rpm.
 - e) Make sure the canal depth is determined by measurement or X-ray verification.
- The planing of the root surface, done with the root facer. Cutting depth: approx. 1 mm, Recommended speed: approx. 2,000 rpm, with little pressure.
- Wind the screw post into the canal with the socket wrench. Move the head-master post back and forth several times while producing the thread.
- a) Clean and disinfect the root canal.
 - b) Fill the canal with cement.
 - c) Cover the post shank with cement.
 - d) Slowly wind in the screw post.
- Shape and prepare the tooth with plastic material such as composites. Finally, build up the restoration.

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GB Pin Master, active Instruction for Use:

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- Treat the root of the tooth endodontically.
 - Remove as much the filling of the root as possible.
 - With the root canal drill prepare the root.
 - Drill the canal to the desired depth. Recommended speed: approx. 500-1000 rpm.
 - Make sure the canal depth is determined by measurement or X-ray verification.
- The root post itself has a self-cutting thread. Use the socket wrench to cut the thread. Move forward slowly, by winding back and forth several times.
- Wind the post out completely and rinse the canal well.
 - Fill the canal with cement. Dip the root post into the cement and insert it into the canal.
 - Remove excess cement with a probe, before it hardens.
- Use Retopin Titan Pins for additional retention.
- Shape and prepare the tooth with plastic material such as composites. Finally, build up the restoration.

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