EDENTA AG, Hauptstrasse 7, CH-9434 AU/SG – Switzerland

Tel.: +41 71 747 25 25 Fax: +41 71 747 25 50 e-mail: info@edenta.ch

EDENTA Etabl., Industriestr. 13, LI-9486 Schaanwald – Liechtenstein

Tel.: +423 375 20 50 Fax: +423 375 20 55 e-mail: info@edenta.com

DE

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GB

Head Master passiv Instruction for Use:

For Direct Restoration With Plastic Filling Materials. System has to be disinfected, cleaned, dried and sterilised before the first use; with steam (autoclave), temperature of 134°C, pressure of 2.1 bar, exposure time 5 min.

Detailed recommendations for preparation in accordance with DIN EN ISO 17664 www.edenta.com

1. a) Treat the root of the tooth endodontically.
b) Remove as much the filling of the root as possible.
c) With the root canal drill prepare the root.
d) Drill the canal to the desired depth. Recommended speed: approx. 500-1000 rpm.
e) Make sure the canal depth is determined by measurement or X-ray verification.
2. The planing of the root surface, done with the root facer. Cutting depth: approx. 1 mm, Recommended speed: approx. 2,000 rpm, with little pressure.
- 3/4. a) As a test, insert the titanium post with the socket wrench, without pressure.
b) Clean and disinfect the root canal.
c) Fill the canal with cement.
d) Cover the post shank with cement.
e) Insert post slowly and carefully into the root canal.
5. Shape and prepare the tooth with plastic material such as composites. Finally, build up the restoration.

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GB

Post Master passiv Instruction for Use:

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Detailed recommendations for preparation in accordance with DIN EN ISO 17664

www.edenta.com

1. a) Treat the root of the tooth endodontically.
b) Remove as much the filling of the root as possible.
c) With the root canal drill prepare the root.
d) Drill the canal to the desired depth. Recommended speed: approx. 500-1000 rpm.
e) Make sure the canal depth is determined by measurement or X-ray verification.
2. a) The cylindrical, passive threaded Titanium post must be tried in the canal.
b) Then the post is shortened to the desired length.
3. a) The post master is cemented.
b) Use Retopin Titan Pins for extra retention.
- 4/5. Shape and prepare the tooth with plastic material such as composites. Finally, build up the restoration.

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