EDENTA AG, Hauptstrasse 7, CH-9434 AU/SG – Switzerland
Tel.: +41 71 747 25 25 Fax: +41 71 747 25 50 e-mail: info@edenta.ch

EDENTA Etabl., Industriestr. 13, LI-9486 Schaanwald – Liechtenstein
Tel.: +423 375 20 50 e-mail: info@edenta.com

DE Anwendungshinweise:

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GB

Instruction for Use:

Retention screw pins for building - up broken - down teeth, for amalgam, composite and crown reconstruction.

1. Select a suitable drilling site in healthy, caries-free dentine, avoiding any damage to the pulp or pericementum. With the drill in the package, drill once until the depth-limit and take it out! Max. 1.000 rpm.
 2. Put the retention screw pin in a right angle and run at approximately 750 RPM. When the pin reaches the drill-hole end, the pin automatically breaks off exactly at the predetermined place. For especially difficult-to-access places, or weak dentine, we recommend inserting the pin manually with the blue handle provided.
 3. If you need more pins for the anchoring of the tooth-filling, repeat the procedure. In some cases it may be necessary to bend the pin for these constructions. The Reto-minipin cannot be bent.
 4. Before the lying of the composite filling, we recommend that the pin should be covered with white or tooth-color.
- The Retopin and Minipin shanks are reusable; simply screw in a new pin from the refill-box. However, the remainder must be removed from the shank by twisting it out.

SAFETY PRECAUTIONS:

- avoidance of pulp vicinity imperative!
- sterilisation PRIOR to use!
- general hygiene instructions to be observed
- consider an approximate distance of 1 mm to the enamel
- use a suitable bending instrument
- apply masking lacquer in case of composite cappings

REMARKS / STORAGE:


In dry conditions and protected against contaminants.
Protect instruments in general against chemicals, acids, heat and extreme temperature variations.

STERILISATION:

autoclave 134°C / 274°F - 2.1 bar - 5 min.

Detailed recommendations for preparation in accordance with DIN EN ISO 17664

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