EDENTA AG, Hauptstrasse 7, CH-9434 AU/SG – Switzerland
Tel.: +41 71 747 25 25 Fax: +41 71 747 25 50 e-mail: info@edenta.ch

EDENTA Etabl., Industriestr. 13, LI-9486 Schaanwald – Liechtenstein
Tel.: +423 375 20 50 e-mail: info@edenta.com

DE Anwendungshinweise Head Master, active:

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GB Head Master, active Instruction for Use:

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- Treat the root of the tooth endodontically.
 - Remove as much the filling of the root as possible.
 - With the root canal drill prepare the root.
 - Drill the canal to the desired depth. Recommended speed: approx. 500-1000 rpm.
 - Make sure the canal depth is determined by measurement or X-ray verification.
- The planing of the root surface, done with the root facer. Cutting depth: approx. 1 mm, Recommended speed: approx. 2,000 rpm, with little pressure.
- Wind the screw post into the canal with the socket wrench. Move the head-master post back and forth several times while producing the thread.
- Clean and disinfect the root canal.
 - Fill the canal with cement.
 - Cover the post shank with cement.
 - Slowly wind in the screw post.
- Shape and prepare the tooth with plastic material such as composites. Finally, build up the restoration.

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GB Pin Master, active Instruction for Use:

For Direct Restoration With Plastic Filling Materials. System has to be disinfected, cleaned, dried and sterilised before the first use; with steam (autoclave), temperature of 134°C, pressure of 2.1 bar, exposure time 5 min.

Detailed recommendations for preparation in accordance with DIN EN ISO 17664 www.edenta.com

1. a) Treat the root of the tooth endodontically.
b) Remove as much of the filling of the root as possible.
c) With the root canal drill prepare the root.
d) Drill the canal to the desired depth. Recommended speed: approx. 500-1000 rpm.
e) Make sure the canal depth is determined by measurement or X-ray verification.
2. a) The root post itself has a self-cutting thread. Use the socket wrench to cut the thread. Move forward slowly, by winding back and forth several times.
b) Wind the post out completely and rinse the canal well.
3. a) Fill the canal with cement. Dip the root post into the cement and insert it into the canal.
b) Remove excess cement with a probe, before it hardens.
4. Use Retopin Titan Pins for additional retention.
5. Shape and prepare the tooth with plastic material such as composites. Finally, build up the restoration.

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