CeraStar & StarTec System

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CeraStar & StarTec System

Premium Tools for Perfect Surfaces on All-Ceramic Materials

Modern all-ceramic restorations place the highest demands on precision, surface quality and material-friendly processing. **The latest generation of CeraStar grinders** enables precise machining across the full spectrum of modern high-performance ceramics – from highly translucent anterior materials to high-strength posterior ceramics with **flexural strengths of 400–1200 MPa**. For controlled, efficient and material-friendly processing at the highest level.

Thanks to the innovative premium hybrid grit and the specially developed hybrid ceramic bond, the CeraStar instruments provide controlled, low-vibration and material-friendly grinding with high cutting efficiency and excellent service life. The combination of diamond and ceramic reduces thermal stress, minimizes microcracks and creates exceptionally homogeneous surfaces – the ideal foundation for subsequent high-gloss polishing.

Optimized for use with CeraStar grinders, high-gloss polishing is performed using the StarTec diamond polishing system. The two-step polishing concept enables fast, safe and material-friendly high-gloss polishing with minimal surface roughness – without any additional glazing cycle.

As an **optimally coordinated system, CeraStar grinders and StarTec polishers** enable efficient machining and high-gloss polishing of all modern all-ceramic materials, including zirconia (ZrO₂), lithium disilicate, zirconia-reinforced lithium silicate (ZLS), hybrid ceramics, CAD/CAM ceramics, veneering ceramics and feldspathic ceramics – delivering highly aesthetic, reproducible and cost-efficient results in both laboratory and clinical applications.

Perfectly coordinated processing system:

1. CeraStar – Precision grinding for all-ceramic materials

CeraStar

opt. 8.000
6.000 - 10.000
max. 5.000*
1

L mm	13,0	11,0	3,0	7,0	2,0	7,0	2,0	8,0	11,0	4,0
Size $\varnothing^{1/10}$ mm	050	040	170	130	130	035	050	040	035	040
Order No.	601.050HP	602.040HP	603.170HP	604.130HP	653.130HP	655.035HP	656.050HP	657.040HP	658.035HP	659.040HP
ISO No. 865 104...	107 512 050	173 512 040	372 512 170	024 512 130	372 512 130	248 512 035	010 512 050	198 512 040	161 512 035	001 512 040

2. StarTec – Efficient high-gloss polishing

StarTec

1

L mm	2,5	14,0	2,0	1,5
Size $\varnothing^{1/10}$ mm	170	040	260	210
opt.	10.000	12.000	5.000	5.000
max.	15.000	15.000	5.000	8.000
Order No.	ST1030HP	ST2030HP	ST1530HP	ST2530HP
ISO No. 804 104...	372 524 170	243 524 040	303 524 260	099 524 210
Stufe / Step / Etape	1			

StarTec

1

L mm	2,5	14,0	2,0	1,5
Size $\varnothing^{1/10}$ mm	170	040	260	210
opt.	7.000	12.000	5.000	5.000
max.	15.000	15.000	5.000	8.000
Order No.	ST1040HP	ST2040HP	ST1540HP	ST2540HP
ISO No. 804 104...	372 514 170	243 514 040	303 514 260	099 514 210
Stufe / Step / Etape	2			

Technical Information


CeraStar – Premium hybrid grit for processing all-ceramic materials

Monolithic restorations – the modern standard in prosthetics

Monolithic restorations combine stability, precision, and durability within a single material structure. They are used in both anterior and posterior regions. Modern all-ceramic materials mainly differ in their flexural strength.

Flexural strength ranges of modern all-ceramic materials

- **Translucent / low-strength ceramics:** 400–700 MPa, ideal for anterior crowns and highly aesthetic applications
- **Medium-strength ceramics:** 600–900 MPa, offering a balanced combination of stability and aesthetics, suitable for anterior restorations and short-span posterior bridges
- **High-strength ceramics:** 900–1200 MPa, ideal for posterior bridges and high-load areas

CeraStar – Universal for all-ceramic materials (Shapes see page 1)

The CeraStar grinders are specially designed for processing the entire spectrum of modern high-performance ceramics – from highly translucent anterior materials to high-strength posterior ceramics with flexural strengths of 400–1200 MPa.

Suitable for:

- Zirconia (ZrO₂)
- Lithium disilicate
- Zirconia-reinforced lithium silicate (ZLS)
- Hybrid ceramics and CAD/CAM ceramics
- Veneering ceramics and feldspar ceramics

Homogeneous surfaces for perfect polishing

CeraStar grinders create uniform, smooth surfaces without irregular grinding grooves, providing the ideal basis for subsequent high-gloss polishing with the StarTec HP polishing system (see page 1).

Benefits:

- Reduced polishing effort
- Uniform material removal
- High-quality, brilliant final surfaces
- Maximum precision and aesthetics

Typical applications

- Shaping and contouring zirconia crowns and bridges
- Precise reduction of lithium disilicate and ZLS restorations
- Fine finishing of veneers, inlays, and overlays
- Chairside adjustments without material damage
- Preparation for subsequent polishing and final finishing

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