


818

Rad
Wheel
Roue



Schaft / Shank	L mm	ISO	Order No.	0,7	0,7	0,7
314 FG standard						
		806 314 042 524...	818.314...	035	040	050
		● 806 314 042 534...	G 818.314...	035	040	050
204 RA standard		806 204 042 524...	818.204...		040	

820


Rad
Wheel
Roue



Schaft / Shank	L mm	ISO	Order No.	2,0	2,0
314 FG standard					
		806 314 043 524...	820.314...	040	050
		● 806 314 043 544...	SG 820.314...	040	050
204 RA standard		806 204 043 524...	820.204...	040	

822


Birne
Pear
Poire



L mm	ISO	Order No.	2,0	2,0	2,0	3,0	
		806 314 237 524...	822.314...	008	009	010	012
		● 806 314 237 514...	F 822.314...				012
		● 806 314 237 534...	G 822.314...	008	009	010	012
		● 806 314 237 544...	SG 822.314...				012

824


Rillenschleifer
Groove Grinder
Instrument de rainurage



L mm	ISO	Order No.	0,8	0,8	
		806 314 055 524...	824.314...	037	042
		Schnitttiefe T mm, cutting depth		0,8	1,0

825


Linse
Knife Edge
Lentille



L mm	ISO	Order No.	0,8	0,8	0,9	0,9	1,5	1,5	
		806 314 304 524...	825.314...	016	018	023	025	040	050
		● 806 314 304 534...	G 825.314...		018	023	025	040	050

830


Birne
Pear
Poire



Schaft / Shank	L mm	ISO	Order No.	2,5	2,7	2,7	2,7	2,7	2,7	2,7	2,7	5,1
314 FG standard												
		806 314 288 524...	830.314...	008	009	010	012	014	016	018	021	
		● 806 314 288 504...	C 830.314...			010	012	014	016			
		● 806 314 288 514...	F 830.314...	008		010	012	014	016			021
		● 806 314 288 534...	G 830.314...			010	012	014	016	018	021	
		● 806 314 288 544...	SG 830.314...				012	014	016	018	021	
313 FG short		806 313 288 524...	S 830.313...				012					

830L

Birne lang
Pear long
Poire, long



Schaft / Shank	L mm	ISO	Order No.	4,0	5,0	5,0	5,0	5,0	7,0
314 FG standard									
		806 314 239 524...	830L.314...	012	014	016	018	021	025
		○ 806 314 239 494...	UF 830L.314...	012					
		● 806 314 239 504...	C 830L.314...		014				
		● 806 314 239 514...	F 830L.314...	012	014	016	018	021	025
		● 806 314 239 534...	G 830L.314...	012	014	016	018	021	025
		● 806 314 239 544...	SG 830L.314...		014	016	018	021	025

834 834A 868A

Tiefenmarkierer
Depth Marker
Marqueur de profondeur



L mm	ISO	Order No.	6,0	6,0	3,5	8,0	8,0	
		806 314 552 524...	834.314...	834	834	834A	868A	868A
		806 314 000 524...	834A.314...	016	021			
		806 314 000 524...	868A.314...			031		
		Schnitttiefe Tmm, cutting depth		0,3	0,5	1,0	0,3	0,4